

Windows 7 – Camera Install

- After extracting the .zip file, the NET2860 Standalone will be in the top directory.
 - Plug camera into USB port
 - Go to “My Computer” – right click
 - Scroll to “Manage” – left click – open device manager
 - Right hand screen scroll to sound video game controllers
 - Right click on USB2860 device update driver software
 - Select: Browse My Computer driver software
 - Select: Let me pick from a list of device drivers on My Computer
 - 32-Bit: Select “Have disk”. Select “Browse” ... NET2860 Standalone\Install_3.0.0.5. Select “NET_2860”. Select “OK”. Select “Next”.
 - 64-Bit: Select “Have disk”. Select “Browse” ... NET2860 Standalone\Install_3.0.0.9_X64. Select “NET_2860_X64”. Select “OK”. Select “Next”. Check the “Always Trust” box, OK.
 - Camera will now be in Device Manager → Imaging device NET_2860 device
 - To install capture: → NET2860 Standalone\Install_1.0.0.0 folder. Run “InstallSnap.exe”. Change Snap Source to “Switch 1”. Click “Install”.
 - You must restart your computer now.
 - Go to NET2860 Standalone and open “VideoView” to capture.
 - We suggest placing a copy of VideoView on your desktop for easy access.
-